# ERGO TOUCH

Multi-Language Instructions For Use  
ErgoTouch FLY / ErgoTouch GLIDE / ECO FLY

Date of Issue: 2021

Version: 2.0

ODC Industries  
Rue Thomas Edison  
ZA Bois Vert  
Ploërmel 56800  
France

CE 0050

## ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE

### Single-use Lens Injection System for injection of One-Piece Foldable Acrylic Lenses.

## Instructions For Use

### RANGE OF MODELS AND COMPONENTS:

Model	Cartridge Size	Body/ Soft	Cartridge	Stopper	Plunger	Spring
ErgoTouch FLY 1.8	1.8mm	SMMA/ Elastomer/ Polypropylene	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch FLY 2.0	2.0mm	SMMA/ Elastomer/ Polypropylene	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch FLY 2.2	2.2mm	SMMA/ Elastomer/ Polypropylene	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch FLY 2.4	2.4mm	SMMA/ Elastomer/ Polypropylene	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch FLY 2.6	2.6mm	SMMA/ Elastomer/ Polypropylene	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch FLY Rail 2.2	2.6mm	SMMA/ Elastomer/ Polypropylene	Polypropylene	Silicone	POM	Stainless Steel
ECO FLY 1.8	1.8mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel
ECO FLY 2.0	2.0mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel
ECO FLY 2.2	2.2mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel
ECO FLY 2.4	2.4mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel
ECO FLY 2.6	2.6mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel
ECO FLY Rail 2.2	2.2mm	SMMA/ Elastomer/ Polypropylene	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch GLIDE 1.8	1.8mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch GLIDE 2.0	2.0mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch GLIDE 2.2	2.2mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch GLIDE 2.4	2.4mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel
ErgoTouch GLIDE 2.6	2.6mm	SMMA	Polypropylene	Silicone	POM	Stainless Steel

### INTEND OF USE

The **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE** injection system facilitates the implantation by an ophthalmic surgeon of a foldable one-piece Hydrophilic and Hydrophobic intraocular lens through incisions of between 1.8 to 2.6 mm, depending on the **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE** model and operating technique used. With the **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE** injection system the majority of all foldable, one-piece acrylic lenses can be injected, provided that the safe injection of the lens using the particular **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE** injection system is validated in advance by the lens manufacturer according to ISO 11979-3 fold and recovery requirements.

The **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE** injection system consists of a sterile single-use injector with silicone cushion and the single-use cartridge.

### INDICATIONS

- Folding and injecting of one-piece foldable lenses into the capsular bag or sulcus after extra-capsular cataract extraction.
- Recommended viscoelastic Sodium ha 1.0 to 1.6 % and compatible with most other viscoelastic concentration (up to 3%). Also compatible with HPMC. The choice is relevant to the ophthalmic surgeon technique
- The choice of cartridge size must be in adequacy with the dioptric power and lens thickness

### CONTRAINDICATIONS

- Three-piece foldable, accommodative and silicone intraocular lenses are **NOT** suited for implantation using the **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE** injection system.
- The **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE** injector must not be used where an IOL manufacturer states a specific injector
- Do not use without Viscoelastic solution.

### INSTRUCTIONS FOR USE

1. Open the blister/pouch pack using standard sterile procedures. Place contents onto a sterile field.
2. Fill the cartridge tunnel of the **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE** cartridge and the loading chamber of the cartridge sufficiently with validated viscoelastic solution. Balanced salt solution should not be used as the sole lubricant.
3. Depending to the used injection system, follow its relevant instructions below:

#### **ErgoTouch FLY / ECO FLY**

- A. Place the lens symmetrically on the central hinge of the open cartridge. Asymmetric insertion of the lens can lead to rotation of the lens in the cartridge.
- B. Press the lens into both guide rails with rounded sterile tweezers and gently move the lens back and forth horizontally in order to ensure that it can move freely.
- C. Ensure that the lens is correctly aligned and make sure that neither the edge of the lens nor the haptic are wedged in.
- D. Press the sides of the cartridge together until the click-lock mechanism engages. The audible click demonstrates the correct closure of the cartridge. A final visual check is recommended after closure to ensure haptic is not trapped.
- E. Insert the cartridge into the front end of the injector. Push the cartridge to the adequate position.

For ECO FLY injector, push the cartridge from the tip of the injector to the very back position until locking (cartridge must be fully locked by the holding parts of the injector).

#### **ErgoTouch GLIDE**

- A. The lens is loaded flat by the back of cartridge. The cartridge has a design to help the tweezer to go in.
- B. The lens should move free of resistance and pushed close to the tunnel cartridge entry.
- C. Ensure that the lens is correctly aligned and make sure that neither the edge of the lens nor the haptic are wedged in.
- D. A final visual check is recommended to ensure the last haptic is pushed close the optic edge.
- E. The cartridge must be inserted by the top of the body. The cartridge design avoids assembling the cartridge on a wrong way into the body nose. A Click must be audible to be sure it is well positioned.

- Carefully push the plunger forward and ensure that the silicone tip correctly enters the loading chamber. If the silicone tip cannot be introduced into the loading chamber, retract the plunger to the starting position and align the silicone tip straight, using sterile tweezers. Continue to push the cushion until the inner spring begins to press together. Pull the plunger back a few millimeters and then push it forward again. This step ensures that the lens is correctly grasped.
- In order to avoid excess viscoelastic entering the eye, gently push the plunger until the lens is at the front of the cartridge, releasing the excess viscoelastic well away from the eye before inserting the injector into the incision.
- Guide the tip of the cartridge through the incision and push it across the iris to the near edge of the pupil.
- Slowly press the plunger forward to support the correct delivery of the lens into the eye and simultaneously withdraw the instrument from the eye. To avoid swelling of the silicone tip during ejection from the cartridge, only push the plunger until the lens has emerged completely.
- Position the lens carefully and if necessary rotate it with the aid of the appropriate positioning hook.
- Remove the viscoelastic material from the eye and the lens, with standard irrigation and aspiration techniques.

#### **OPERATIONAL PROCEDURE**

The appropriate surgical techniques and eligibility of the patients are the responsibility of the ophthalmic surgeon. He or she must assess the appropriateness of the relevant procedure based on his or her experience, and the directions for use with the **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE**.

#### **PRECAUTIONS**

- The **ErgoTouch FLY / ECO FLY / ErgoTouch GLIDE** should be used by trained and qualified professionals only. The result of surgical procedure depends on skill and abilities used with the device and other products during the surgery.
- This product is supplied sterile and should not be used if inner packaging is soaked, opened or damaged.
- All contents are for single use and should never be re-sterilized.
- Package must be opened according to required aseptic technique.
- Do not use the cartridge if tip has been broken or damaged as this may causes scratches or damage to the Lens.
- Always use viscoelastic when loading the foldable intraocular lens into the cartridge.
- Please follow up carefully the loading instructions.
- In case the device falls on the floor or touch any other unsterile product, discard it according to standard rules for residues in health establishments.

If the above precautions are not followed, foreseeable risks or complications could cause:

- Improper handling, loading and use of the injector can cause damage to the lens and patient such as tissue damage, inflammation, contamination
- Contamination during use could generate infection / contamination of the patient
- Re-use of the Single-use Injector and/or cartridge could generate contamination between patients

#### **GUARANTEE AND LIABILITY LIMITATION**

The manufacturer guarantees that this product has been manufactured with appropriate care and assumes no liability for side effects or resulting damages, losses or costs that may arise as a result of the direct or indirect use of this product. Manufacturer's liability is restricted to the performance of repairs resulting from product defects, which are clearly not the result of incorrect handling or the use of lenses not validated with this injector model.

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